**OPED GmbH**  
Medizinpark 1  
83626 Valley/Oberlaindern  
Germany  
Fon +49 (0) 80 24/60 81 8-210  
Fax +49 (0) 80 24/60 81 8-299  
mail@oped.de  
www.oped.de

**OPED AG**  
Hinterbergstrasse 26  
6330 Cham  
Switzerland  
Fon +41 (0) 41/748 53 35  
Fax +41 (0) 41/748 53 40  
contact@oped.ch  
www.oped.ch

**OPED UK Ltd**  
Unit 2 Glenmore Business Centre  
Waller Road  
Hopton Industrial Estate  
Devizes  
SN10 2EQ  
United Kingdom  
Phone +44 (0) 1380 722177  
Fax +44 (0) 1380 710095  
enquiries@oped.biz  
www.opeduk.co.uk

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Ihr OPED-Team

## OPTIVOshoulder

Intended use: orthosis to immobilise the shoulder joint in an adjustable position.

### Safety Advice

 The product is single use and should only be used by one patient. 

Should it be used for the treatment of more than one patient, the manufacturer's product liability shall lapse in accordance with the Medical Devices Act.

 Never apply the product to open wounds and only use under prior medical instruction. 

 The patient may only adjust it in consultation with his/her treating doctor.

 Should you experience excessive pain, swelling, numbness or sensitivity changes, unusual reactions or an uncomfortable sensation while

wearing the SUPROshoulder, please contact your doctor or orthopaedic technician immediately.

 This orthosis was developed for the indications stated and may be damaged through unusual and non-intended use or through activities such as sport, falls or other accidents, as well as liquids.

 Please consult your treating doctor before using it in conjunction with another medical device.

 Please observe the road transport law of your country.

 Note: Not made with natural rubber latex. 

### How to Use

#### Application with 15° abduction pads

● Fasten abdominal belt (a). Insert the arm support (b) into the abduction pad (c) until you hear a snap. Move sideways towards the OPTIVOshoulder and place the arm in the arm support (b). (Fig. 1)

● Fasten velcro strap. Lay shoulder strap (d) over the healthy shoulder. (Fig. 2-3) Attach hook (e) and position the belt clip (f). (Fig. 4-7) The length of the shoulder strap can be adjusted using the green belt clasp (k) in the shoulder pad. (Fig. 8)

#### Application with 0° abduction

● Fasten abdominal belt (a). Insert the hook pad (g) with the green strap into the ringlet (h) and the hook (i) into the arm support. (Fig. 9-11) Move sideways towards the OPTIVOshoulder and place the arm in the arm support (b). (Fig. 12)

● Fasten velcro strap. Lay shoulder strap (d) over the healthy shoulder. (Fig. 13-14) Attach hook (e) and position the belt clip (f). (Fig. 15-16) Attach the OPTIVOshoulder to the abdominal belt in the desired position using the hook pad (g). (Fig. 17)

### Arm support length adjustment

- Press the button down and adjust to the desired setting (j).

### Hand support length adjustment

- Press the clips on the bottom together and slide them into the desired position. The hand support is detachable (k).

### Cleaning and Care

**Washing:** wash strap and pad in the washing machine at 40°C

Clean plastic parts with a damp cloth and with plastic-friendly cleaning agents only.

**Drying:** leave to dry naturally

### Product Characteristics

Sizes: Sizes: S/M and L/XL

### Indications

- Dislocations
- Humeral head fracture
- Clavicle fracture
- Scapula fracture
- AC joint separation
- Rupture of the rotator cuff
- Supraspinatus tendon
- Endoprostheses

### Contraindications

None known

### Technical data/parameters

- Material composition:  
Arm support: PP, small parts: PA, POM  
strap, pad: PP, PS, cotton.
- Dimensions:  
(without strap): 44cm x 22cm x 21cm
- Weight: 1050 g

### Transport conditions/ storage conditions/disposal:

**Storage:** Store in a dry place and out of direct sunlight.

### Guarantee

The guarantee is in line with statutory prescriptions.

### Patient information

**Showering/hygiene:**  
You can shower/be in water while wearing the product, subject to the approval of your treating doctor.

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